

# 標準液槽 ChromTech® Standard Cells

## 標準石英液槽 (兩面磨光,兩面磨砂)

適用波長：190nm-2500nm

產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	每組(2支)售價
 標準石英液槽(含蓋)	45 x 12.5 x 3.5	1	10	0.35	Q11001	
	45 x 12.5 x 4.5	2	10	0.70	Q11002	
	45 x 12.5 x 7.5	5	10	1.7	Q11005	
	<b>45 x 12.5 x 12.5</b>	<b>10</b>	<b>10</b>	<b>3.5</b>	<b>Q11010</b>	
	45 x 12.5 x 22.5	20	10	7.0	Q11020	
	45 x 12.5 x 32.5	30	10	10.5	Q11030	
	45 x 12.5 x 42.5	40	10	14	Q11040	
	45 x 12.5 x 52.5	50	10	17.5	Q11050	
 標準石英液槽(含塞)	45 x 12.5 x 7.5	5	10	1.7	Q12005	
	<b>45 x 12.5 x 12.5</b>	<b>10</b>	<b>10</b>	<b>3.5</b>	<b>Q12010</b>	
	45 x 12.5 x 22.5	20	10	7.0	Q12020	
	45 x 12.5 x 32.5	30	10	10.5	Q12030	
	45 x 12.5 x 42.5	40	10	14.0	Q12040	
	45 x 12.5 x 52.5	50	10	17.5	Q12050	
	45 x 12.5 x 102.5	100	10	35.0	Q12100	

## 標準玻璃液槽 (兩面磨光,兩面磨砂)

適用波長：340nm-2500nm



產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	每組(4支)售價
 標準玻璃液槽(含蓋)	45 x 12.5 x 3.5	1	10	0.35	G11001	
	45 x 12.5 x 4.5	2	10	0.70	G11002	
	45 x 12.5 x 7.5	5	10	1.7	G11005	
	<b>45 x 12.5 x 12.5</b>	<b>10</b>	<b>10</b>	<b>3.5</b>	<b>G11010</b>	
	45 x 12.5 x 22.5	20	10	7.0	G11020	
	45 x 12.5 x 32.5	30	10	10.5	G11030	
	45 x 12.5 x 42.5	40	10	14	G11040	
	45 x 12.5 x 52.5	50	10	17.5	G11050	
 標準玻璃液槽(含塞)	45 x 12.5 x 7.5	5	10	1.7	G12005	
	<b>45 x 12.5 x 12.5</b>	<b>10</b>	<b>10</b>	<b>3.5</b>	<b>G12010</b>	
	45 x 12.5 x 22.5	20	10	7.0	G12020	
	45 x 12.5 x 32.5	30	10	10.5	G12030	
	45 x 12.5 x 42.5	40	10	14.0	G12040	
	45 x 12.5 x 52.5	50	10	17.5	G12050	
	45 x 12.5 x 102.5	100	10	35.0	G12100	

# 螢光液槽 ChromTech® Fluorometer Cells

## 螢光石英液槽 (四面磨光)



產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	每組(2支)售價
 螢光石英液槽(含蓋)	45 x 12.5 x 7.5	5	10	1.7	Q21005	
	<b>45 x 12.5 x 12.5</b>	<b>10</b>	<b>10</b>	<b>3.5</b>	<b>Q21010</b>	
	45 x 12.5 x 22.5	20	10	7.0	Q21020	
	45 x 12.5 x 32.5	30	10	10.5	Q21030	
	45 x 12.5 x 42.5	40	10	14	Q21040	
	45 x 12.5 x 52.5	50	10	17.5	Q21050	
	45 x 12.5 x 7.5	5	10	1.7	Q21005	
 螢光石英液槽(含塞)	45 x 12.5 x 7.5	5	10	1.7	Q22005	
	<b>45 x 12.5 x 12.5</b>	<b>10</b>	<b>10</b>	<b>3.5</b>	<b>Q22010</b>	
	45 x 12.5 x 22.5	20	10	7.0	Q22020	
	45 x 12.5 x 32.5	30	10	10.5	Q22030	
	45 x 12.5 x 42.5	40	10	14.0	Q22040	
	45 x 12.5 x 52.5	50	10	17.5	Q22050	
	45 x 12.5 x 7.5	5	10	1.7	Q22005	

## 螢光微量石英液槽 (四面磨光)

產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	每組(2支)售價
 螢光微量石英液槽(含蓋)	<b>45 x 12.5 x 12.5</b>	<b>10</b>	<b>2</b>	<b>0.7</b>	<b>Q23102</b>	
	<b>45 x 12.5 x 12.5</b>	<b>10</b>	<b>4</b>	<b>1.4</b>	<b>Q23104</b>	
 螢光微量石英液槽(含塞)	<b>45 x 12.5 x 12.5</b>	<b>10</b>	<b>2</b>	<b>0.7</b>	<b>Q24102</b>	
	<b>45 x 12.5 x 12.5</b>	<b>10</b>	<b>4</b>	<b>1.4</b>	<b>Q24104</b>	

# 微量液槽 ChromTech® Micro Cells

## 微量石英液槽 (黑邊避光)

產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	每組(2支)售價
 黑邊微量石英液槽(含蓋)	45 x 12.5 x 12.5	10	1	0.35	Q31101	
	45 x 12.5 x 12.5	10	2	0.7	Q31102	
	45 x 12.5 x 12.5	10	3	1	Q31103	
	45 x 12.5 x 12.5	10	4	1.4	Q31104	
	45 x 12.5 x 12.5	10	5	1.75	Q31105	
 黑邊微量石英液槽(含塞)	45 x 12.5 x 12.5	10	2	0.7	Q32102	
	45 x 12.5 x 12.5	10	4	1.4	Q32104	

## 微量石英液槽 (磨砂白邊)

產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	每組(2支)售價
 白邊微量石英液槽(含蓋)	45 x 12.5 x 12.5	10	1	0.35	Q41101	
	45 x 12.5 x 12.5	10	2	0.70	Q41102	
	45 x 12.5 x 12.5	10	3	1	Q41103	
	45 x 12.5 x 12.5	10	4	1.4	Q41104	
	45 x 12.5 x 12.5	10	5	1.7	Q41105	
 白邊微量石英液槽(含塞)	45 x 12.5 x 12.5	10	2	0.7	Q42102	
	45 x 12.5 x 12.5	10	4	1.4	Q42104	

# 避光流動液槽 ChromTech® *Self Masking Flow*

## 流動石英液槽 (黑邊避光)

產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	每組(2支)售價
 避光流動石英液槽	45 x 12.5 x 12.5	10	φ1.5	0.018	Q51105	
	45 x 12.5 x 12.5	10	φ2	0.03	Q51102	
	45 x 12.5 x 12.5	10	φ3	0.07	Q51103	
 避光流動石英液槽	45 x 12.5 x 7.5	5	4x12	0.24	Q52054	
	45 x 12.5 x 12.5	10	4x12	0.48	Q52104	
	45 x 12.5 x 22.5	20	4x12	1.0	Q52204	
	45 x 12.5 x 32.5	30	4x12	1.5	Q52304	
 避光流動石英液槽	45 x 12.5 x 12.5	10	φ1	0.01	Q53101	
	45 x 12.5 x 12.5	10	φ1.5	0.018	Q53105	
	45 x 12.5 x 12.5	10	φ2	0.03	Q53102	

## 其他液槽



Micro cell with black walls



Cylindrical cell



Large cell



Micro cell with frosted walls